

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# PCDGF,GRN Fusion Protein



## Basic Information

**Catalog Number:**

Ag6198

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MWTLVSWVALTAGLVAGTRCPDGGQFCPVACCLDPG  
GASYSCCRPLLDKWPTTISRHLGGPCQVDAHCSAGH  
SCIFTVSGTSSCCPFPEAVACGDGHHCCPRGFHCSA  
DGRSCFQRSGNNSVGAIQCPDSQFECPDFSTCCVM  
VDGSWGCCPMPQASCCEDRVHCCPHGAFCDLVHT  
RCITPTGTHPLAKKLPARTNRAVALSSVMCPDARS  
RCPDGSTCCELPSGKYGCCPMPNATCCSDHLHCCPQ  
DTVCDLIQSKLSKENATDILLTKLPAHTVGDVKCDM  
EVSCPDPGYTCCRLQSGAWGCCPFTQAVCCEDHIHC  
CPAGFTCDTQKGTCEQPHQVPWMEKAPAHLSLPD  
PQALKRDVPCDNVSSCPSSDTCRDRNRQGWACCPY  
RQGVCCADRRHCCPAGFRCAARGTKCLRREAPRWD  
APLRDPALRQLL

(1-593aa encoded by BC000324)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.