

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# BIN1 Fusion Protein

## Basic Information

### Catalog Number:

Ag6240

### Form:

Available lyophilized

### Species:

human

### Expression Source:

*e. coli*-derived, PGEX-4T, with N-terminal GST.

### Biological Activity:

Not tested

### Endotoxin Level:

Please contact the lab for more information

### Validated Application:

Blocking peptide

### Peptide Sequence:

MAEMGSKGVTAGKIASNVQKILTRAQEKVLQKLGKA  
DETKDEQFEQCVQNFNKQLTEGTRLQKDLRTYLASV  
KAMHEASKKLNELQEVYEPDWPGRDEANKIAENND  
LLWMDYHQKLVLDQALLTMDTYLGQFPDIKSRIAKRGR  
KLVVDYDSARHHYESLQTAKKKDEAKIAKAEELIKAQK  
VFEEMNVDLQEELPSLWNSRVGFYVNTFQSIAGLEEN  
FHKEMSKLNQNLNDVLVGLEKQHGSNFTVKAQPSD  
NAPAKGNKSPSPDGSPAATPEIRVNHEPEPAGGATP  
GATLPKSPSQPAEASEVAGGTQPAAGAQPGETAAS  
EAASSSLPAVVVETFPATVNGTVEGGSGAGRLDPPG  
FMFKVQAQHDYATDDELQLRAGDVVLVLPFNPE  
EQDEQWLMGVKESDWNQHKLEKCRGVFPENFTE  
VP

(1-439 aa encoded by BC004101)

## Reconstitution and Storage

### Reconstitution:

Reconstitute at 0.25 µg/ µl in 200 µl sterile water for short-term storage.

After reconstitution with sterile water, if glycerol has no effect on subsequent experiments, it is recommended to add an equal volume of glycerol for long-term storage (see Stability and Storage for more details).

If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used).

Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.

### Shipping:

The product is shipped at ambient temperature.

Upon receipt, store it immediately at the recommended temperature (see below).

## Stability and Storage

Store for up to 12 months at -20°C to -80°C as lyophilized powder.

## Storage of Reconstituted Protein


### Short Term Storage:

Store at 2-8°C for (1-2 weeks).

### Long Term Storage:

Aliquot and store at -20°C to -80°C for up to 3 months, reconstitution with sterile water and addition of an equal volume of glycerol. Avoid repeat freeze-thaw cycles.

## Selected Validation Data



For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.