

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SEC14L2 Fusion Protein



Basic Information

Catalog Number:

Ag6267

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSGRVGDLSPRQKEALAKFRENVQDVLPALPNPDDY
FLLRWLRARSFDLQKSEAMLRKHVEFRKQKDIDNIIISW
QPPEVIQQYLSGGMCGYDLGGCPVWYDIIGPLDAKG
LLFSASKQDLLRTKMRECELLQECAHQTTKLGRKVET
ITIIYDCEGLGLKHLWKPAVEAYGEFLCMFEENYPETLK
RLFVVKAPKLPVAYNLIKPFLESDTRKKIMVLGANWK
EVLLKHISPDQVPVEYGGTMDPDGPNKCKSKINYGG
DIPRKYVVRDQVKQQYEHSVQISRGSSHQVEYEILFP
GCVLRWQFMSDGDVGFGLFKTKMGERQRAGEMT
EVLPNQRYNSHLVPEDGTLTCSDPGICKYLC LGNALKP
HVQLSACEVPLPPWIFGSEC
(1-392aa encoded by BC058915)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.