

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant mouse Elmod1 protein



## Basic Information

Catalog Number:

Ag6325

Size:

50 µg

Form:

Available lyophilized

Species:

mouse

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MKHFLRMLIQVCLYFYCKFLWRCKMFKVMRKLTGRCEL  
QRICVGTKPGASRTMKIETSLRDSKSKLLQTSVSVHPD  
AIEKTIIDDIMELKKINPDINPQLGISLQACLLQIVGYRNL  
IADVEKLRREPYSDNPQHEEMLLKLWELLKPNTPLES  
RVSKQWCEIGFQGDDPKTDFRGMGLLGLYNLQYFAE  
RDATVAQQVLSDSVHPKCSKFSKIEWEKKKMDKAIG  
YSFAIVGINITDLAYNLLVSGALKTHFYNIAPAPTLSH  
FQQTFCYLMHEFHKFWIEEDPMDIMEFNRVREKFRKR  
IIKQLQNPDMALCPHFAASEGLINM  
(1 - 326 aa encoded by NM\_177769)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.