

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

CNPY2,MSAP Fusion Protein



Basic Information

Catalog Number:

Ag6363

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MKGWGWLALLGALLGTAWARRSQDLHCGACRALV
DELEWEIAQVDPKKTIQMGFSFRINPDGSQSVVEVPY
ARSEAHLTELEEICDRMKEYGEQIDPSTHRKNYVRV
GRNGESSELDQGI RIDS D I S G T L K F A C E S I V E E Y E D E L I
E F F S R E A D N V K D K L C S K R T D L C D H A L H I S H D E L
(1-182 aa encoded by BC065015)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.