

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ERG protein



## Basic Information

Catalog Number:

Ag6402

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

QWLEWAVKEYGLPDVNILFQNI DGKELCKMTKDDF  
QRLTPSYNADILLSHLHYLRETPLPHLTSDDV D KALQN  
SPRLMHARNTGGAAFI F PNTSVYPEATQRITTRPDL PY  
EPPRRSAWTGHGHPTQSKAAQ P SPSTVPKTEDQRP  
QLDPYQILGPTSSRLANPGSGQIQLWQF LLELLSDSSN  
SSCITWEGTNGEFKMTDPDEVARRWGERKSKPNMNY  
DKLSRALRYYYDKNIMTKVHGKRYA YKDFH G I A QAL  
QPHPPESLYKYPDLPMGSYHAHPQKMN F V A P H P  
PALPVTSSSFFAAPNPY WNSPTGGIYPNTRLPTSHMP  
SHLGTYY

(141-479 aa encoded by BC040168)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.