

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant mouse Fig4 protein



Basic Information

Catalog Number:

Ag6429

Size:

50 µg

Form:

Available lyophilized

Species:

mouse

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

KMQSCYEKME~~SL~~RLDGLQQRFDVSS~~TS~~VFKQRAQIH
MREQMDNAVYTFETLLHQELGKGP~~TK~~EELCKSIQRVL
ERV~~L~~KKYDYDSSSV~~R~~KRFFREALLQISIPFLLK~~L~~APTCK
SELPRFQELIFEDFARFILVENTYEEV~~L~~LQTVMKDILQ
AVKEAAVQRKHNLYRDSMVMHNSDPN~~L~~HLLAEGAPI
DWGEEYSNSGGGGSPSP~~ST~~PESATLSEKRRRAKQVVS
VVQDEEVGLPFEASPE~~SP~~PPASPDGVTEIRGLLAQGL
RPESPPAGPLLNGALAGESPQPKAAPEASSPPASPL
QHLLPGKAVDLGPPKPSDQETGEQVSSPSSHPALHTT
TEDSAGVQTE

(111 - 456 aa encoded by BC015295)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.