

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human UBOX5 protein



Basic Information

Catalog Number:

Ag6573

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

QPAKTCSEQVIDSILLVTSENLPQDVALQAPALPMES
DCDPGDQPESQAPSSLQKLAELIQDVPEEFDPITLE
IMPCPMLLP SGKVIDQSTLEKCNREATWGRVPSDPF
TGVAFTPHSQPLPHPSLKARIDHFLQHSIPGCHLLGR
AQ TALAVIPSSIVLP SQKRKIEQAEHVPDSNFGVNAS
CFSATSPVLPTTSEHTAKKMKATNEPSLTHMDCSTG
PLSHEQKLSQSLEIALASTLGSMPSFTARLTRGQLQHL
GTRGSNTSWRPGTGSEQPGSILGPECASCKRVFSPYF
KKEPVYQLPCGHLLCRPCLGEKQRSLPMTCTACQRPV
ASQDVLRVHF

(196-541 aa encoded by BC000515)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.