

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# SPCS2 Fusion Protein



## Basic Information

Catalog Number:

Ag6640

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MAAAVQGGRRSGGSGGCSGAGGASNCGTGSGRSG  
LLDKWKIDDKPVKIDKWDGSAVKNSLDDSVKVLLEK  
YKYVENFGLIDGRLTICTISCFFAIVALIWDYMHPPES  
KPVLAALCVISYFVMMGILTIYTSYKEKISFLVAHRKDPT  
GMDPDDIWQLSSSLKRFDDKYTLKLTIFSGRTKQQR  
AEFTKSIKFFDHSGLVMDAYEPEISRLHDSLAIERKI  
K

(1-226 aa encoded by BC064957)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.