

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

NAP1L1 Fusion Protein



Basic Information

Catalog Number:

Ag6690

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MADIDNKEQSELDQDLDDVEVEVEEETGEETKLKARQ
LTVQMMQNPQILAALQERLDGLVETPTGYIESLPRVV
KRRVNALKNLQVKCAQIEAKFYEEVHDLERKYAVLYQ
PLFDKRFEIIINAIYEPTEECEWKPDEEDEISEELKEKAK
IEDEKDEEKEDPKGPEFWLTVFKNVLLSDMVQEH
DEPILKHLKDIKVKFSDAGQPMSFVLEHFEPNEYFTN
EVLTKTYMRSEPDSDPF5FDGPEIMGCTGCQIDWK
KGNVTLKTIKKKQKHKGRGTVRTVTKTVSNDFFNF
FAPPEVPESGDLDDDAEAILAADFEIGHFLRERIIPRSV
LYFTGEAIEDDDDDYDEEGEEADEEGEEEGDEENDPD
YDPKKDQNPAECKQQ

(1-391 aa encoded by BC002387)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.