

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human TMC6 protein



Basic Information

Catalog Number:

Ag6767

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, T-HIS, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MIQSPEAAGQEAVLLVLPVVGLLNLGAPYLCRVLAAL
EPHDSPLVLEVYVAICRNLIKLAILGTLCYHWLGRRVG
VLQGGQWEDFVGQELYRFLVMDVFLMLLDTLFGELV
WRIISEKCLKRRRKPEFDIARNVLELYGQTLTWLGVLF
SPLLPAVQIIKLLLVFYVKKTSLLANCQAPRRPWLASH
MSTVFLTLFCFPAFLGAAVFLCYAVWQVKPSSTCGPF
RTLDTMYEAGRWWRHLEAAGPRVSWLPWWHRYLM
ENTFFVFLVSALLLAVIYLNIVVRGQRKVICLLKEQIS
NEGEDIKFLINKLHSIYERKEREERSRVGTTEEAAPPA
LLTDEQDA

(C-term-345aa encoded by BC018346)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.