

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SYK Fusion Protein



Basic Information

Catalog Number:

Ag6773

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6[°]His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

RKLLTLEDKELGSGNFGTVKKGYQMKKVVKTVAVKIL
KNEANDPALKDELLAEANVMQQLDNPYIVRMIGICEA
ESWMLVMEMAELGPLNKYLQNRHVKDKNIIELVHQ
VSMGMKYLEESNFVHRDLAARNVLLVTQHYAKISDFG
LSKALRADENYYKAQTHGKWPVKWYAPECINYYKFSS
KSDVVSFGVLMWEAFSYGQKPYRGMKGSEVTAMLE
KGERMGCPAGCPREMYDLMNLCWTYDVENRPGFAA
VELRLRNYYYDVVN

(367-635 aa encoded by BC001645)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.