

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human DHODH protein



## Basic Information

Catalog Number:

Ag6810

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MAWRHLKRAQDAVILGGGGLLFASYLMATGDERFY  
AEHLMPTLQGLLDPESAHRNAVRFSLGLLPRARFQDS  
DMLEVRVLGHKFRNPVGIAGFDKHGEAVDGLYKMG  
FGFVEIGSVTPKPKQEGNPRPRVFRLPEDQAVINRYGF  
NSHGLSVVEHLRARQKQAKLTEDGLPLGVNLGKNK  
TSVDAAEYAEGVRLGPLADYLVVNVSSPNTAGLRS  
LQGKAE LRRLTKVLQERDGLRRVHRPAVLVKIAPDLT  
SQDKEDIASVVKELGIDGLIVTNTTVSRPAGLQ GALRS  
ETGGLSGKPLRDLSTQTIREMYALTQGRVPIIGVGGVS  
SGQDALEKIRAGASLVQLYTALTFWGPVVGKVKREL  
EALLKEQGFGGVTDAIGADHRR

(1-395 aa encoded by BC065245)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.