

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human AKNA protein



Basic Information

Catalog Number:

Ag6912

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

RLED^{SL}HQPLQGSPT^{RP}ASAFDRPARTRGRPADSPAT
WGS^{HY}GSKSTERLPGEP^{RG}EEQIVPPGRQ^{RA}RSSSVP
REVLRL^{SL}SESELPSL^{PL}SEKSK^{TT}KD^{SP}QAARDGKR
GVGSAGWPDRV^{TR}FGQYTGHE^{YH}VLSPKAVPKGNG
TVSCPHCRPI^{RT}QDAGGAV^{TD}PLGPPADTLQCPL
CGQVGS^{PE}ADGPGSATSGAEKAT^{TR}RRKAPST^{PS}PK
QRSKQAGSS^{PR}PPGLWYLATAPPAP^{PA}FAYISSV
PI^{MP}YP^{PA}AV^{YY}APAG^{PT}SAQ^{PA}AK^{WP}PTAS^{PP}PAR
RHRHSIQ^{LD}LG^{DL}EELNKALSRAV^{QA}AESV^{RS}TRQM
RSSLSADLRQA^{HS}LRG^SCLF

(1090-1439aa encoded by BC055285)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.