

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RASSF9 Fusion Protein



Basic Information

Catalog Number:

Ag6919

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

SERVLPLTRILKLWKAWGDEQPNMQFVLVKADAFLP
VPLWRTAEAKLVQNTKELWELSPANYMKTLPDPKQKR
IVRKTRKLAIKIKQDTVSHDRDNMETLVHLIISQDHTI
HQQVKRMKELDLEIEKCEAKFHLDREVNDGENYVQD
AYLMPSESEVEQNLDLQYEENQTLLEDLSESDGIEQLEE
RLKYYRILIDKLSAEIEKEVKSVCIDINEDAEGEAAASELE
SSNLESVKCDLEKSMKAGLKIHSHLSGIQKEIKYSDSLL
QMKAKEYELLAKEFNLSLHISNKDGCQLKENRAKESEV
PSSNGEIPPFQRFVFSNYTNDTSDTGISSNHSQDSE
TTVGDVLLST

(C-term-350aa encoded by BC031589)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.