

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# HNRNPF Fusion Protein



## Basic Information

**Catalog Number:**

Ag7128

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MMLGPEGGEGFVVKLRGLPWSCSVEDVQNFLSDCTI  
HDGAAGVHFIYTRGRQSGEAFVELGSEDDVKMALK  
KDRESMGHRYIEVFKSHRTEMWV LKHSGPNSADSA  
NDGFVRLRGLPFGCTKEEIVQFFSGLEIVPNGITLPVD  
PEGKITGEAFVQFASQELAEKALGKHKERIGHRYIEVF  
KSSQEEVRSYSDPPLKFMSVQRPGPYDRPGTARRYIGI  
VKQAGLERMRPGAYSTGYGGYEEYSGLSDGYGFTTDL  
FGRDLSYCLSGMYDHRYG DSEFTVQSTTGHCVHMARG  
LPYKATENDIYNFFSPLNPRVHIEIGPDGRVTGEADV  
EFATHEEAVAAMSKDRANMQHRYIELFLNSTTGASN  
GAYSSQVMQGMGVSAAQATYSGLESQSVSGCYGA  
GYSQQNSMGGYD

(1-415 aa encoded by BC001432)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.