

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human GALK1 protein



## Basic Information

**Catalog Number:**

Ag7300

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

MAALRQPQVAELLAERRAFREEFGAEPELAVSAPGR  
VNLIGEHTDYNQGLVLPMALELMTVLVGSPRKDGLVS  
LLTSEGADEPQRLQFPLPTAQRSLEPGTPRWANYVK  
GVIQYYPAAPLPGFSAVVVSSVPLGGGLSSASLEVA  
TYTFLQQLCPDSGTIAARAQVCQQAHSFAGMPCGI  
MDQFISLMGQKGHALLIDCRSLETSLVPLSDPKLAVLIT  
NSNVRHSLASSEYPVRRRQCEEVARALGKESLREVQL  
EELEAARDLVSKEGFRRARHVVGIEIRTAQAAAALRR  
GDYRAFGRMLVESHRSRDDYEVSCPELDQLVEAALA  
VPGVYGSRMTGGGFGGCTVTLEASAAPHAMRHIQE  
HYGGTATFYLSQAADGAKVLCL

(1-392 aa encoded by BC001166)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.