

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

EFCAB4B Fusion Protein



Basic Information

Catalog Number:

Ag7358

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MAAPDGRVSRPQRLGQSGQGPKGSGACLHPLDSL
EQKETQEQTSGQLVMLRKAQEFFQTCDAEGKGFAR
KDMQRLHKELPLSLEELEDVFDALDADGNGYLTPOEF
TTGFSHFFSQNNPSQEDAGEQVAQRHEEKVYLSRG
DEDLGDMGEDEEAQFRMLMDRLGAQKVLEDES DVK
QLWLQLKKEEPHLLSNFEDFLTRISQLQEAHEEKNELE
CALKRKIAAYDEEIQHLYEEMEQQIKSEKEQFLLKDTE
RFQARSQEQELKLLCKEQELEQLTQKQKRLGQCTAL
HHDKHETKAENTKLLTNQELARELERTSWELQDAQ
QQLESLQEQACKLHQEKEMEVYRVTESLQREKAGLLK
QLDFLRCVGGHWPVLRAPPSLGSEGPV
(1-395 aa encoded by BC004524)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.