

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human B4GALT4 protein



Basic Information

Catalog Number:

Ag7517

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

GAIQEI PKAKEFMANFHKTLLGKGTLTNEASTKKVEL
DNCPSVSPYLRGQSKLIFKPDLTLEEVAENPKVSRGR
YRPEECKALQRVAIILVPHRNREKHLMYLLEHLHPFLQR
QQLDYGIYVIHQAEKKFNRAKLLNVGYLEALKEENW
DCFIHVDLVPENDFNLYKCEEHPKHLVGRNSTGY
RLRYSGYFGGVTALSREQFFKVNNGFSNNYWGWGGED
DDLRLRVELQRMKISRPLPEVGKYTMVFHTRDKGNEV
NAERMKLLHQVSRVWRTDGLSSCSYKLVSEHNPLYI
NITVDFWFGA

(36-344 aa encoded by BC004523)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.