

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human UBAP2L protein



Basic Information

Catalog Number:

Ag7715

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

TSSTVSAPPPVSVSSSLNSGSSGLSLGSNSTVTAST
RSSVATTSGKAPPNLP PGVPPLLPNPYIMAPGLLHAYP
PQVYGYDDLQMLQTRFPLDYYSIPFPTPTTPLTGRDGS
LASNPYSGDLTKFGRGDASSPAPATTLAQPPQNNQTQ
THHTTQQTFNLNPPGYSYTSLPYYTGVPGLPSTFQY
GPAVFPVAPTSSKQHGVNVS VNASATPFQQPSGYG
SHGYNTGVSVTSSNTGVDPISGSVYSKTQQSFEKQG
FHSGTPAASFNLPSALGSGGPINPATAAAYPPAPFMH
ILTPHQQPHSQILHHHLQDGGQTSGGQRSQTSSIPQ
KPQTNKSAYNYSYSGAN

(739-1087 aa encoded by BC003170)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.