

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human GRINA protein



Basic Information

Catalog Number:

Ag7785

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MSHEKSFLVSGDNYPNPGYGGPQPPMPPYAQPP
YPGAPYPQPPFQSPYQPGYPHGSPYPQGGYPQG
PYPQGGYPQGPYPQEGYPQGPYPQGGYPQGPYPQN
PFPPNPYGPQVFPQDDPDSPOHGNYQEEGPPSY
DNQDFPATNWDDKSIRQAFIRKVFVLTQLSVTLST
VSVFTFVAEVKGFVRENVWTYYVSYAVFFISLIVLSC
CGDFRRKHPWNLVALSVLTASLSYMGMIASFYNTEA
VIMAVGITTAVCFTVIFSMQTRYDFTSCMGVLLVS
MVLVFIFAILCIFIRNRILEIVYASLGALLFTCF LAVDTQ
LLGNKQLSLSPEEYVFAALNLYTDIINIFLYLTIIGRAK
E

(1-371aa encoded by BC041788)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.