

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

WDR5 Fusion Protein



Basic Information

Catalog Number:

Ag7900

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MA~~T~~EEK~~P~~PETE~~A~~ARAQ~~P~~TPSS~~S~~ATQ~~S~~KTP~~V~~KPN~~Y~~AL
K~~F~~TLAG~~H~~T~~K~~AV~~S~~SV~~K~~FSP~~N~~GE~~W~~LASS~~S~~AD~~K~~L~~I~~K~~I~~W~~G~~A
Y~~D~~G~~K~~F~~E~~K~~T~~I~~S~~G~~H~~K~~L~~G~~I~~S~~D~~V~~A~~W~~S~~S~~D~~S~~N~~L~~L~~V~~S~~A~~S~~D~~D~~K~~T~~L~~K~~
I~~W~~D~~V~~S~~S~~G~~K~~C~~L~~K~~T~~L~~K~~G~~H~~S~~N~~Y~~V~~F~~C~~C~~N~~F~~N~~P~~Q~~S~~N~~L~~I~~V~~S~~G~~S~~F
D~~E~~S~~V~~R~~I~~W~~D~~V~~K~~T~~G~~K~~C~~L~~K~~T~~L~~P~~A~~H~~S~~D~~P~~V~~S~~A~~V~~H~~F~~N~~R~~D~~G~~S~~L~~I
V~~S~~S~~S~~Y~~D~~G~~L~~C~~R~~I~~W~~D~~T~~A~~S~~G~~Q~~C~~L~~K~~T~~L~~I~~D~~D~~D~~N~~P~~P~~V~~S~~F~~V~~K~~F~~S
P~~N~~G~~K~~Y~~I~~L~~A~~A~~T~~D~~N~~T~~L~~K~~L~~W~~D~~Y~~S~~K~~G~~K~~C~~L~~K~~T~~Y~~T~~G~~H~~K~~N~~E~~K~~Y~~C
I~~F~~A~~N~~F~~S~~V~~T~~G~~G~~K~~W~~I~~V~~S~~G~~S~~E~~D~~N~~L~~V~~Y~~I~~W~~N~~L~~Q~~T~~K~~E~~I~~V~~Q~~K~~L~~Q
G~~H~~T~~D~~V~~V~~I~~S~~T~~A~~C~~H~~P~~T~~E~~N~~I~~I~~A~~S~~A~~A~~E~~N~~D~~K~~I~~K~~L~~W~~K~~S~~D~~C~~
(1-334 aa encoded by BC001635)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.