

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ADAMTS18 protein



Basic Information

Catalog Number:
Ag8038

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
KRPA Y T W S I V Q S E C S V S C G G G Y I N V K A I C L R D Q N T Q
V N S S F C S A K T K P V T E P K I C N A F S C P A Y W M P G E W S T C
S K A C A G G Q Q S R K I Q C V Q K K P F Q K E E A V L H S L C P V S T
P T Q V Q A C N S H A C P P Q W S L G P W S Q C S K T C G R G V R K
R E L L C K G S A A E T L P E S Q C T S L P R P E L Q E G C V L G R C P K N
S R L Q W V A S S W S E C S A T C G L G V R K R E M K C S E K G F Q G
K L I T F P E R R C R N I K K P N L D L E E T C N R R A C P A H P V Y N M V
A G W Y S L P W Q Q C T V T C G G G V Q T R S V H C V Q Q G R P S S
S C L L H Q K P P V L R A C N T N F C P A P E K R E D P S C V D F F N W
C H L V P Q H G V C N H K F Y G K Q C C K S C T R K I
(872-1221aa encoded by BC063283)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.