

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human BYSL protein



## Basic Information

**Catalog Number:**  
Ag8052

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Peptide Sequence:**  
MPKFKAARGVGGQEKHAPLADQILAGNAVRAQVREK  
RRGRGTGEAEEYVGPRLSRRILQARQQEELEAEH  
GTGDKPAAPRETRRLGPRMPQDGSDEDEEWPTLE  
KAATMTAAGHHAEEVVDPEDERAIEFMNKNPPAR  
RTLADII MEKLTEKQTEVETVMSEVSGFPMPQLDPRVL  
EYYRGVREVLSKYRSGIKPKAFKIIIPALSNWEQILYVTE  
PEAWTAAAMYQATRIFASNLKERMAQRFYNLVLLPR  
VRDDVAEYKRLNFHLYMALKKALFKPGA WFKGILPLC  
ESGTCTLREAIIVGSIITKCSIPVLHSSAAMLKIAEMEYS  
GANSIFLRLLLDKKYALPYRVL DALVFHFLGFRTEKRELP  
VLWHQC LLTLVQRYKADLATDQKEALLELLRLQHPHQ  
(1-437aa encoded by BC050645)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.