

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# BYSL Fusion Protein



## Basic Information

**Catalog Number:**

Ag8052

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MPKFKARGVGGQEKHAPLADQILAGNAVRAGVREK  
RRGRGTGEAEEYVGPRLSRRLQARQQEELEAEH  
GTGDKPAAPRETRRLGPRMPQDGSDEDEEWPTLE  
KAATMTAAGHHAEEVVDPEDERAIEFMNKNPPAR  
RTLADII MEKLTEKQTEVETVMSEVSGFPMPQLDPRVL  
EYYRQVREVL SKYRSGKLPKAFKIIPALSNWEQILYVTE  
PEAWTAAAMYQATRIFASNLKERMAQRFYNLVLLPR  
VRDDVAEYKRLNFHLYMALKKALFKPGA WFKGILPLC  
ESGTCTLREAIIVGSIITKCSIPVLHSSAAMLKIAEMEYS  
GANSIFLRLLLDKKYALPYRVL DALVFHFLGFRTKRELP  
VLWHQC LLTLVQRYKADLATDQKEALLELLRLQPHPQ  
(1-437aa encoded by BC050645)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.