

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# RTN4IP1 Fusion Protein



## Basic Information

Catalog Number:

Ag8100

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MEFLKTCVLRNACTAVCFWRSKVVQKPSVRRISTTS  
PRSTVMPAWVIDKYGKNEVLRFTQNMMPHIIHYPNE  
VIVKVHAASVNPIDVNMRSYGATALNMKRDP LHVKI  
KGEEFPLTLGRDVSQVMECGLDVKYFKPGDEVWAA  
VPPWKQGTLSEFVVVSGNEVSHKPKSLTHTQAASLP  
YVALTAWSAINKVGLNDKNCCTGKRVLILGASGGVGT  
FAIQVMKAWDAHVAVCSQDASELVRKLGADDVID  
YKSGSVEEQLKSLKPFDFILDVGGSTETWAPDFLKK  
WSGATYVTLVTPFLNMDRLGIADGMLQTVTVGSK  
ALKHFVWGVHYRWAFFMASGCPCLDDIAELVDAGKIR  
PVIEQTFPFSKVPEAFLKVERGHARGKTVINNV  
(1-396aa encoded by BC006399)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.