

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# BAT4 Fusion Protein



## Basic Information

**Catalog Number:**

Ag8350

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MSRPLLI FT PAT DP SD L WK D G Q Q Q P Q P E K P E S T L D G  
A A A R A F Y E A L I G D E S S A P D S Q R S Q T E P A R E R K R K K R R I  
M K A P A A E A V A E G A S G R H G Q G R S L E A E D K M T H R I L R A  
A Q E G D L P E L R R L L P H E A G G A G G N I N A R D A F W W T P L  
M C A A R A G Q G A A V S Y L L G R G A A W V G V C E L S G R D A A  
Q L A E E A G F P E V A R M V R E S H G E T R S P E N R S P T P S L Q Y C  
E N C D T H F Q D S N H R T S T A H L L S L S Q G P Q P P N L P L G V P I  
S S P G F K L L R G G W E P G M G L G P R G E G R A N P I P T V L K R D  
Q E G L G Y R S A P Q P R V T H F P A W D T R A V A G R E R P P R V A T  
L S W R E E R R R E E K D R A W E R D L R T Y M N L E F  
(1-356 aa encoded by BC008783)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.