

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human AGXT2L2 protein



Basic Information

Catalog Number:

Ag8374

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MGKSI¹GN²GH³PV⁴AC⁵VAA⁶TQP⁷VAR⁸A⁹FEAT¹⁰GV¹¹EY¹²FNT¹³
GG¹⁴SPV¹⁵SCAV¹⁶GLAV¹⁷LN¹⁸VLE¹⁹KE²⁰QL²¹QD²²HAT²³SV²⁴GS²⁵FL²⁶M²⁷Q²⁸L²⁹
LG³⁰QQ³¹KIK³²HP³³IV³⁴GD³⁵VR³⁶GV³⁷GL³⁸FIG³⁹VD⁴⁰L⁴¹IK⁴²DE⁴³AT⁴⁴RT⁴⁵P⁴⁶A⁴⁷E⁴⁸E⁴⁹
AA⁵⁰Y⁵¹LV⁵²S⁵³R⁵⁴L⁵⁵KEN⁵⁶Y⁵⁷LL⁵⁸ST⁵⁹D⁶⁰GP⁶¹GR⁶²N⁶³IL⁶⁴K⁶⁵F⁶⁶K⁶⁷PP⁶⁸MC⁶⁹F⁷⁰SL⁷¹D⁷²NA⁷³
R⁷⁴Q⁷⁵V⁷⁶V⁷⁷AK⁷⁸LD⁷⁹AIL⁸⁰TD⁸¹M⁸²E⁸³E⁸⁴K⁸⁵V⁸⁶R⁸⁷SC⁸⁸ET⁸⁹LR⁹⁰L⁹¹Q⁹²
(276-450 aa encoded by BC008009)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.