

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

AGXT2L2 Fusion Protein



Basic Information

Catalog Number:

Ag8374

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MGKSI¹GN²GH³PV⁴AC⁵VA⁶AT⁷QP⁸VAR⁹A¹⁰FE¹¹AT¹²GV¹³EY¹⁴F¹⁵N¹⁶T¹⁷
GG¹⁸SP¹⁹V²⁰SC²¹AV²²GL²³AV²⁴LN²⁵VLE²⁶KE²⁷QL²⁸QD²⁹HAT³⁰SV³¹GS³²FL³³M³⁴Q³⁵L³⁶
LG³⁷QQ³⁸K³⁹IK⁴⁰HP⁴¹IV⁴²GD⁴³VR⁴⁴GV⁴⁵GL⁴⁶FI⁴⁷GV⁴⁸DL⁴⁹IK⁵⁰DE⁵¹AT⁵²RT⁵³PA⁵⁴TE⁵⁵E⁵⁶
AA⁵⁷YL⁵⁸VS⁵⁹R⁶⁰L⁶¹KEN⁶²YV⁶³LL⁶⁴ST⁶⁵DG⁶⁶PGR⁶⁷N⁶⁸IL⁶⁹K⁷⁰FP⁷¹PM⁷²CF⁷³SL⁷⁴DN⁷⁵A⁷⁶
RQ⁷⁷VV⁷⁸AK⁷⁹LD⁸⁰AIL⁸¹TD⁸²MEE⁸³KV⁸⁴RSC⁸⁵ET⁸⁶LR⁸⁷L⁸⁸Q⁸⁹P⁹⁰
(276-450 aa encoded by BC008009)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.