

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ANT3 protein



## Basic Information

Catalog Number:

Ag8521

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MTEQAI<sup>S</sup>FAKDFLAGGIAAAISKTAVAPIERVKLLQV  
QHASKQIAADKQYKGI<sup>V</sup>DCIVRIPKEQGVLSFWRGNL  
ANVIRYFPTQALNFAFKDKYKQIFLGGVDKHTQFWRY  
FAGNLASGGAAGATSLCFVYPLDFARTRLAADV<sup>G</sup>KS  
GTEREFRGLGDCLVKITKSDGIRGLYQGFVSVQGI<sup>I</sup>Y  
RAAYFGVYDTAKGMLPDPK<sup>N</sup>THIVVSWMIAQTVTAV  
AGVVSYPFD<sup>T</sup>VRRRMMMQSGRKGADIMYGT<sup>V</sup>DCW  
RKIFRDEGGKAFFKGAWSN<sup>V</sup>LRGMGGAFVLLVLYDELK  
KVI

(1-298aa encoded by BC008935)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.