

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human FGL1 protein



Basic Information

Catalog Number:

Ag8593

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MAKVFSFILVTTALTMGREISALEDCAQEQMRLRAQV
RLETRVKQQQVKIKQLQENEVQFLDKGDENTVIDL
GSKRQYADCSEIFNDGYKLSGFYKIKPLQSPAEFVYC
DMSDGGGWTVIQRRSDGSENFNRGWKDYENGFNGF
VQKHGEYWLGNKNLHFLTTEQEDYTLKIDLADFEKNSR
YAQYKNFKVGVDEKFNFYELNIGEYSGTAGDSLGNFHP
EVQWWASHQRMKFSTWDRDHDNYEGNCAEEDQS
GWWFNRCHSANLNGVYYSGPYTAKTDNGIVWYTWLH
GWWYSLKSVVMKIRPNDFIPNVI
(1-312 aa encoded by BC007047)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.