

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# BHMT Fusion Protein



## Basic Information

**Catalog Number:**

Ag8738

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MPPVGGKKAKKILERLNAGEIVIGDGGFVFALEKRGY  
VKAGPWTPEAAVEHPEAVRQLHREFLRAGSNVMQT  
FTFYASEDKLENRGNVLEKISGQEVNEAACDIARQV  
ADEGDALVAGGVSQTPSYLSCKSETEVKKVFLOQLE  
VFMKKNVDFLIAEYFEHVVEAVWAVETLIASGKPVAA  
TMCIGPEGDLHGVPPGECVRLVKAGASII GVNCHFD  
PTISLKTVKLMKEGLEAARLKAHLSQPLA YHTPDCNK  
QGFIDLPEFPF GLEPRVATRWDIQKYAREAYNLGVRYI  
GGCCGFEPYHIRAIAEELAPERGF LPPASEKHGWSGS  
GLDMHTKPWVRARARKEYWENLRIASGRPYNPSMSK  
PDGWGVTKGTAEMLMQQKEATTEQQLKELFEKQKFKS  
Q

(1-406 aa encoded by BC012616)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.