

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# LIPT1 Fusion Protein



## Basic Information

Catalog Number:

Ag8825

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6<sup>His</sup>.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

VVIGRHQNPWQECNLNLMREEGIKLARRRSGGGTVY  
HDMGNINLTFFTTKKYYDRMENKLVRLNAVQPQL  
DVQATKRFDLLDGGQFKISGTASKIGRTTAYHHCITLLC  
STDGTFLLSLLKSPYQGI RSNATASIPSLVKNLLEKDPTL  
TCEVLMNAVATEYAA YHQIDNHHILINPTDETLFPGIN  
SKAKELQ TWEWIYGKTPKFSINTSFHVLVEQSHLEIKV  
FIDIKNGRIEICNIEAPDHWLPLEIRDKLNSSLIGSKFCP  
TETMLTNILLRTPQDHLKNSKWNILCEKIKGIM  
(72-373 aa encoded by BC007001)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.