

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

NUP54 Fusion Protein



Basic Information

Catalog Number:

Ag8882

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MAFNFGAPSGTSGTAAATAAPAGGFGGFGTTSTTA
GSAF5FSAPTNTGTTGLFGGTQNKGF5GFGTGT
GTSTGLGTGLTGLGFGGNTQQQQTTLGGF5SQ
TQAPTQSNQLINTASALSAPTLGDERDAI LAKWNQL
QAFWGTGKGYFNNNIPPVEFTQENPFCRFKAVGYSC
MPSNKDEDGLVVLVFNKKETEIRSQQQLVESLHKVL
GGNQTLTVNVEGKTLPDDQTEVVIYVVERSPNGTSR
RVPATTL5AHFEQANIKTQLQQLGVTLSMTRTELSPA
QIKQLLQNP5AGVDPIIWEQAKVDNPDSEKLI5P5PMV
GFKELLRRLK5VQDQMTKQH5TRL
(1-351 aa encoded by BC012559)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.