

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SERPINB2 Fusion Protein



Basic Information

Catalog Number:

Ag8906

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MEDLCVANTLFA~~N~~LFLKHLAKASPTQNLFLSPW~~S~~ISST
MAMVYMGSRG~~S~~TE~~D~~QMAKVLQFNEVGANAVTPMT
PENFTSCGF~~M~~Q~~Q~~I~~Q~~KGSYPDAI~~L~~QAQAADKIHSSFR
SLSSAINASTGNYLLESV~~N~~KLFGK~~S~~ASF~~R~~E~~E~~YIRLCQK
YYSSEPQAVDFLECAEEARKKIYSWVITQTKGKIPNLL
PEGSVDG~~D~~TRMVLVNAVYFKGKWKTPFEK~~K~~LNGLYPF
RVNSAQRTPVQMMYLREKLNIGYIEDLKAQILELPYAG
DVSMFLLLPDEIADVSTGLELLESEITYDKLNKWT~~S~~KDK
MAEDEV~~E~~VYIPQFKLEE~~H~~YELRSILRSMGMEDAFNKG
RANFSGMSERNDLFLSEV~~F~~HQAMVDVNEEGTEAAA
GTGGVMTGRTGHGGPQFVADHPFLIMHKITNCILF
FGRFSSP

(1-415 aa encoded by BC012609)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.