

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# MAP3K6 Fusion Protein



## Basic Information

Catalog Number:

Ag9311

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e. coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MNLLSYRDVQDYSIIELVETLQALPTCDVAEQHNV  
CFHYTFALNRRNRPGDRAKALSPLLPLVQLEGSVAPDL  
YCMCGRIYKDMFFSSGFQDAGHREQAYHWYRKAFD  
VEPSLHSGINAAVLLIAAGQHFEDSKELRLIGMKLGCLL  
ARKGCVEKMQYYWVDFYLGAIQLANDPTQVVLA  
EQLYKLNAPIWYLVSVMETFLLYQHFRPTPEPPGGPPR  
RAHFWLHFLQSCQPFKTACAQGDQCLLVLEMNK  
VLLPAKLEVRGTDVSTVTLSEPETQDIPSSWTFPV  
ASICGVSASKRD

(1-307 aa encoded by BC015914)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.