

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human SERPINB2 protein



Basic Information

Catalog Number:

Ag9319

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MEDLCVANTLFA~~NL~~FKHLAKASPTQNLFLSPW~~SIS~~ST
MAMVYMGSRG~~STED~~QMAKVLQFNEVGANAVTPMT
PENFTSCGF~~MQ~~QIQKGSYPDAI~~LQAQA~~ADKIHSSFR
SLSSAINASTGNYLLESV~~NKLF~~GEKSASF~~REEY~~IRLCQK
YYSSEPQAVDFLECAEEARKKIYSWV~~IKTQ~~TGKIPNLL
PEG~~SVD~~GDTRMVLVNAVYFKGKWKTPFEK~~KL~~NGLYPF
RVNSAQRTPVQMMYLREKLNIGYIEDLKAQILELPYAG
DVSMFLLLPDEIADVSTGLELLESEITYDKLNK~~WTS~~KDK
MAEDEV~~EVI~~PQFKLEEHYELRSILRSMGMEDAFNKG
RANFSGMSERNDLFLSEV~~FHQ~~AMVDVNEEGTEAAA
GTGGVMTGRTGHGGPQFVADHPFL~~LI~~MHKITNCILF
FGRFSSP

(1-415 aa encoded by BC012609)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.