

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human Alpha-1-Antitrypsin protein



Basic Information

Catalog Number:

Ag9369

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

QGDAAQKTDTSHHDDHPTFNKITPNLAFAFSLYR
QLAHQSNSTNIFFSPVSIATAFAMLSLGTKADTHDEIL
EGLNFNLTIEPAQIHEGFQELLRTLNQPDSQLQLTTG
NGLFLSEGLKLVDFLEDVKKLYHSEFTVNFGDTEEA
KKQINDYVEKGTQGGKIVDLVKELDRDTVFALVNYIFFK
GKWERPFVEVKDTEEDFHVDQATTVKVPMMKRLGM
FNIQHCKKLSSWLLMKYLGNAIIFLPDEGKLQHLE
NELTHDIITKFLNEDRRSASLHLPKLSITGTYDLKSVLG
QLGITKVFSGADLSGVTEEAPLKLSKAVHKAVLTIDE
KGTEAAGAMFLEAIPMSIPPEVKFNKPFVFLMIEQNTK
SPLFMGKVVNPTQK

(28-418 aa encoded by BC015642)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.