

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RNF31 Fusion Protein



Basic Information

Catalog Number:

Ag9440

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, T-HIS, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

RQHFTIALKEKHITDMVCPACGRPDLTDDTQLLSYFST
LDIQLRESLEPDAYALFHKKLTEGVLMRDPKFLWCAQ
CSFGFIYEREQLEATCPQCHQTFVCRCKRQWEEQHR
GRSCEDFQNWKRMDPEYQAQGLAMYQENGIDCP
KCKFSYALARGGCMHFHCTQCRHQFCSGCYNAFYA
KNKCPEPNCRVKKSLSHGHPDCLFYLRDWTALRLQK
LLQDNNVMFNTEPPAGARAVPGGGCRVIEQKEVPNG
LRDEACGKETPAGYAGLCQAHYKEYLVSLINAHSLDP
ATLYEVEELETATERYLHVRPQPLAGEDPPAYQARLLQ
KLTTEEVLGQSIPIRRRK

(258-603 aa encoded by BC017376)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.