

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human BANP protein



Basic Information

Catalog Number:

Ag9480

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

VRC AV P G R R Q N T I V V K V P G Q E D S H H E D G E S G S E A S D
S V S S C G Q A G S Q S I G S N V T L I T L N S E E D Y P N G T W L G D E
N N P E M R V R C A I I P S D M L H I S T N C R T A E K M A L T L L D Y L F
H R E V Q A V S N L S G Q G K H G K K Q L D P L T I Y G I R C H L F Y K F
G I T E S D W Y R I K Q S I D S K C R T A W R R K Q R G Q S L A V K S F S
R R T P N S S S Y C P S E P M M S T P P P A S E L P Q P Q P Q P Q A L H Y
A L A N A Q Q V Q I H Q I G E D G Q V Q V I P Q G H L H I A Q V P Q G
E Q V Q I T Q D S E G N L Q I H H V G Q D G Q V L Q G A Q L I A V A S
S D P A A A G V D G S P L Q G S D I Q V Q Y V Q L A P V S D H T A G A
Q T A E A L Q P T L Q P E M Q L E H G A I Q I Q
(119-469 aa encoded by BC009424)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.