

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

CFI Fusion Protein



Basic Information

Catalog Number:

Ag9905

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MKLLHVFLLFLCFHLRFCKVITYTSQEDLVEIKKCLAKKY
THLSCDKVFCQPWQRCIEGTCVCKLPYQCPKNGTAV
CATNRRSFPTYCQQKSLECLHPGTFLNNGTCTAEGK
FSVSLKHGNTDSEGIVEVKLVQDKTFMICKSSWSMR
EANVACLDLGFQQGADTQRRFKLSDSLINSTECLHVH
CRGLESLAECTFTKRRTMGYQDFADVVCYTQKADS
PMDDFFQCVNGKYISQMKACDGINDCGDQSDLECC
KACQGGKFHCKSGVCIPSYQCNGEVDCITGEDEV
GCAGFASVAQEETEILTADMDAERRRIKSLPKLSCGV
KNRMHIRRKRIVGGKRAQLGKMKQISLDFKGLHRVYA
ILEFCCIIL

(1-376aa encoded by BC020718)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.