

For Research Use Only

# Human IgG Polyclonal antibody

Catalog Number: B900640

1 Publications



## Basic Information

Catalog Number:

B900640

Size:

1 mg/ml

Source:

N/A

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

155 kDa

Purification Method:

Protein A

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

This product is Normal human IgG which purified with Protein A. The concentration is 1.8ug/ul. An isotype control antibody should have the same immunoglobulin type and be used at the same concentration as the test antibody.

## Notable Publications

Author	Pubmed ID	Journal	Application
Junlan Zhong	35200359	Biosensors (Basel)	

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

