

For Research Use Only

Goat anti-Human IgG serum

Catalog Number: B900720



Physical State	Liquid
Source	Collected from Goat which has been immunized by human IgG. Aseptically collected and processed to minimize hemolytic activity.
Size	1ml
Preservative	0.02% Sodium Azide
Storage	-20°C, avoid repeated thawing and freezing
Note	This product is for research use only, not for diagnostic or therapeutic use.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.