

For Research Use Only

Biotin Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody

Catalog Number: Biotin-65571



Basic Information

Catalog Number:

Biotin-65571

Concentration:

100ug, 500 ug/ml

Source:

Mouse

Isotype:

IgG2a, kappa

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

MOPC-173

Recommended Dilutions:

FC: 0.25 ug per 10⁶ cells in a 100 µl suspension

Applications

Tested Applications:

FC

Species Specificity:

n/a

Positive Controls:

FC : human PBMCs,

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 human PBMCs were surface stained with 0.25 μ g Biotin Mouse IgG2a Isotype Control (MOPC-173) RecAb (Biotin-65571, Clone: MOPC-173) and CoraLite®488-conjugated streptavidin (red), or Biotin Anti-Human CD8 (UCHT4) Mouse IgG2a Recombinant Antibody (Biotin-65559, Clone: UCHT4) and CoraLite®488-conjugated streptavidin, or unstained (blue). Cells were not fixed.