

For Research Use Only

# Biotin Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody

Catalog Number: Biotin-65607



## Basic Information

Catalog Number:

Biotin-65607

Size:

100ug, 500 µg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

24796

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

W3/13

## Applications

Tested Applications:

FC

Species Specificity:

rat

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> rat splenocytes were surface stained with 0.25 ug Biotin Anti-Rat CD43 (W3/13) Mouse IgG2a RecAb (Biotin-65607, Clone: W3/13) and CoraLite®488-conjugated streptavidin, or unstained. Cells were not fixed.



1x10<sup>6</sup> rat splenocytes were surface stained with 0.25 ug Biotin Anti-Rat CD43 (W3/13) Mouse IgG2a RecAb (Biotin-65607, Clone: W3/13) and CoraLite®488-conjugated streptavidin (red), or unstained (blue). Cells were not fixed.