

For Research Use Only

Biotin-conjugated LSP1 Recombinant monoclonal antibody

Catalog Number: Biotin-87378



Basic Information

Catalog Number: Biotin-87378	GenBank Accession Number: NM_002339.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 4046	CloneNo.: 252179A7
Isotype: IgG	UNIPROT ID: P33241-1	Recommended Dilutions: IHC: 1:50-1:500
Immunogen Catalog Number: EG5464	Full Name: lymphocyte-specific protein 1	
	Calculated MW: 37 kDa	
	Observed MW: 52 kDa	

Applications

Tested Applications: IHC	Positive Controls: IHC : human tonsillitis tissue,
Species Specificity: human	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Lymphocyte-specific protein 1(LSP1), also known as leukocyte-specific protein 1, is an F-actin binding and cytoskeleton associated protein. The protein is expressed in lymphocytes, neutrophils, and macrophages. LSP1 is a substrate for mitogen-activated protein kinase (MAPK)-activated protein kinase 2 and protein kinase C. The implication of LSP1 in diverse biological processes, such as inflammation and malignancy, suggests the major regulatory role of LSP1 in signal transduction(PMID: 38254711).

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using Biotin-87378 (LSP1 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using Biotin-87378 (LSP1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).