

For Research Use Only

# FcZero-rAb® Biotin Anti-Mouse CD117/c-Kit Rabbit Recombinant Antibody

Catalog Number: Biotin-FcA98041



## Basic Information

Catalog Number:

Biotin-FcA98041

Concentration:

100ug, 500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0883

GenBank Accession Number:

NM\_001122733.1

GeneID (NCBI):

16590

UNIPROT ID:

P05532-1

Full Name:

kit oncogene

Calculated MW:

109 kDa

Purification Method:

Protein A purification

CloneNo.:

240246C11

Recommended Dilutions:

FC: 0.25 ug per 10<sup>6</sup> cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

-

## Applications

Tested Applications:

FC, ELISA

Species Specificity:

mouse

Positive Controls:

FC : mouse bone marrow cells,

## Background Information

CD117, also known as c-Kit and SCFR, is a transmembrane protein with tyrosine kinase activity encoded by the oncogene c-kit (PMID: 2448137). It is a member of the type III receptor tyrosine kinase family, which also includes CSF-1R, PDGFR  $\beta$ , PDGFR  $\alpha$ , and FLT3 (PMID: 29518044). CD117 is expressed on hematopoietic stem cells and progenitor cells, mast cells, and is also found in a wide range of non-haemopoietic cell types (including melanocytes, germ cells, astrocytes, renal tubules, breast glandular epithelial cells, sweat glands, and interstitial cells of Cajal) (PMID: 10582338; 23073628). CD117 plays an important role in early haemopoiesis. It is also involved in pigmentation, fertility, gut movement, and some aspects of the nervous system (PMID: 23073628).

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> mouse bone marrow cells were surface stained with APC Anti-Mouse CD45R, and 0.25 ug Biotin Anti-Mouse CD117/c-Kit Rabbit RecAb (Biotin-FcA98041, Clone: 240246C11) or FcZero-rAb™ Biotin Rabbit IgG Isotype Control Recombinant Antibody (Biotin-FcA98136, Clone: 240953C9), and Streptavidin-PE Conjugate (PE-PF00030). Cells were not fixed.



1x10<sup>6</sup> mouse bone marrow cells were surface stained with 0.25 ug Biotin Anti-Mouse CD117/c-Kit Rabbit RecAb (Biotin-FcA98041, Clone: 240246C11) (red) or FcZero-rAb™ Biotin Rabbit IgG Isotype Control Recombinant Antibody (Biotin-FcA98136, Clone: 240953C9) (blue), and Streptavidin-PE Conjugate (PE-PF00030). Cells were not fixed.