

For Research Use Only

FcZero-rAb® Biotin Anti-Human E-cadherin Rabbit Recombinant Antibody

Catalog Number: Biotin-FcA98123



Basic Information

Catalog Number:

Biotin-FcA98123

Concentration:

100ug, 500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1112

GenBank Accession Number:

NM_004360.5

GeneID (NCBI):

999

UNIPROT ID:

P12830-1

Full Name:

cadherin 1, type 1, E-cadherin (epithelial)

Calculated MW:

97 kDa

Purification Method:

Protein A purification

CloneNo.:

241098F7

Recommended Dilutions:

FC: 0.25 ug per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

-

Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : HT-29 cells,

Background Information

Cadherins are a family of transmembrane glycoproteins that mediate calcium-dependent cell-cell adhesion and play an important role in the maintenance of normal tissue architecture. E-cadherin (epithelial cadherin), also known as CD324, CDH1 (cadherin 1) or CAM 120/80, is a classical member of the cadherin superfamily which also include N-, P-, R-, and B-cadherins. E-cadherin is expressed on the cell surface in most epithelial tissues. The extracellular region of E-cadherin establishes calcium-dependent homophilic trans binding, providing specific interaction with adjacent cells, while the cytoplasmic domain is connected to the actin cytoskeleton through the interaction with p120-, α -, β -, and γ -catenin (plakoglobin). E-cadherin is important in the maintenance of the epithelial integrity, and is involved in mechanisms regulating proliferation, differentiation, and survival of epithelial cell. E-cadherin may also play a role in tumorigenesis. It is considered to be an invasion suppressor protein and its loss is an indicator of high tumor aggressiveness.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ HT-29 cells (left) or Jurkat cells (right) were surface stained with 0.25 ug Biotin Anti-Human E-cadherin Rabbit RecAb (Biotin-FcA98123, Clone: 241098F7) (red) or 0.25ug Isotype Control (blue), and Streptavidin-PE Conjugate (PE-PF00030). Cells were not fixed.