

For Research Use Only

# FcZero-rAb® Biotin Anti-Human CCL22/MDC Rabbit Recombinant Antibody

Catalog Number: Biotin-FcA98280



## Basic Information

Catalog Number:

Biotin-FcA98280

Concentration:

100ug, 500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1974

GenBank Accession Number:

BC027952

GeneID (NCBI):

6367

UNIPROT ID:

O00626

Full Name:

chemokine (C-C motif) ligand 22

Calculated MW:

93 aa, 11 kDa

Purification Method:

Protein A purification

CloneNo.:

241668D5

Recommended Dilutions:

FC (Intra): 0.25 ug per 10<sup>6</sup> cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

-

## Applications

Tested Applications:

FC (Intra), ELISA

Species Specificity:

human

Positive Controls:

FC (Intra): human monocyte-derived mature dendritic cells,

## Background Information

CCL22, also known as macrophage-derived chemokine (MDC), is a member of the CC chemokine family. It is primarily produced by myeloid cells such as macrophages and dendritic cells (DCs) under steady-state conditions. CCL22 plays a significant role in immune responses by regulating the migration of different subsets of leukocytes. It shares the same receptor, CCR4, which is predominantly expressed by Th2 cells, making CCL22 a key mediator in the development of Th2-dominant diseases such as atopic dermatitis and asthma. Elevated serum levels of CCL22 are observed in patients with these conditions. CCL22 is also expressed in several types of tumor cells and Foxp3+ regulatory T cells (Tregs), suggesting its involvement in attracting Tregs into the tumor microenvironment to evade immune attacks.

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> human monocyte-derived mature dendritic cells were intracellularly stained with 0.25 ug Biotin Anti-Human CCL22/MDC Rabbit RecAb (Biotin-FcA98280, Clone: 241668D5) (red) or FcZero-rAb™ Biotin Rabbit IgG Isotype Control Recombinant Antibody (Biotin-FcA98136, Clone: 240953C9) (blue), and Streptavidin-PE Conjugate (PE-PF00030). Cells were fixed and permeabilized with Intracellular Flow Cytometry Fixation & Permeabilization Buffer Kit (PF00019).