

For Research Use Only

Conjugation CleanAb Kit for Rabbit IgG, Agarose



www.ptgcn.com

Catalog Number: CK001

Basic Information

Catalog Number:
CK001

Applications:
Purification

Conjugate:
Agarose beads; bead size: ~ 90 µm (cross-linked 4
% agarose beads)

Host:
/

Type:
/

Class:
Recombinant - Animal free production

Description

Binding capacity /

Elution buffer /

Wash buffer compatibility /

Affinity (K_D) /

Storage

Storage:

Shipped at ambient temperature. Upon receipt store at +4°C. Stable for one year. DO not freeze!

Storage Buffer:
20% Ethanol

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com
W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.